

SciPette

Liquid Handling Solutions

FLOW Pipette Controller



www.sciquip.co.uk

SciPette FLOW

- Comfortable and simple to use
- Precise Pipetting control
- Robust and lightweight
- Fits the entire range of volumetric and graduated pipettes both glass and plastic, from 0.1 to 100ml
- Uses replaceable filter
- Easy to maintain and clean
- Four colours

Pipette Controllers

Part no.	Colour
SP2205-BLU	blue
SP2205-RED	red
SP2205-GRE	green
SP2205-YEL	yellow



Compress the suction bellows before attaching the pipette

Expel residual liquid when using blow-out pipettes

Filling and delivery are easily controlled with the small sensitive lever

Replace the filter and clean the adapter if necessary after unscrewing the adapter support

SciQuip Ltd
Newtown
Wem
Shrewsbury
Shropshire
SY4 5NU, UK

Tel: +44 (0)1939 234222
Fax: +44 (0)1939 234221
Email: info@sciquip.co.uk
www.sciquip.co.uk